

CNAP OEM Module

Making continuous, non-invasive, arterial blood pressure measurement available to all medical disciplines

The CNAP OEM module is designed to expand medical devices of interested manufacturers with CNAP monitoring. With an integrated CNAP OEM module, such devices are capable to measure blood pressure truly continuously (beat-to-beat) and to display and record data usually only available from invasive measurement:

- beat-to-beat systolic, diastolic and mean blood pressure,
- a high fidelity blood pressure curve and
- the pulse rate

A customized version of the CNAP OEM Module is already part of [the Dräger® Infinity® patient monitoring system](#).

For measurement ranges and the specification of the module contact CNSystems' CNAP product specialists. [E-mail](#)

