

Features

CNAP *continuous non-invasive arterial pressure*

- high-fidelity blood pressure waveform
- real time systolic, diastolic and mean blood pressure
- absolute accuracy and autozero using NBP
- ECG high resolution 3-channel:
- 6-channel display synchronized with patient signals (e.g. blood pressure)

Impedance Cardiography - continuous non-invasive cardiac output

- real time beat-to-beat cardiac output and stroke volume
- complete hemodynamic profile including total peripheral resistance

ANS *Autonomic Nervous System*

- autoregulation of blood pressure
- heart rate and blood pressure variability
- baroreceptorreflex sensitivity

Easy-to-use software packages:

- Syncope
- Pacemaker
- Print report / export of data



Task Force® Monitor,
tilted patient & physician